








drea™ Model specifications

	Diameter	Min. Tube Size	Tip angle	Tip length	Cutter type	Shield	Label color code
3S45FU	3mm	10mm	45°	12mm	Fluted	Downward, Partial	
3S45FUS	3mm	10mm	45°	12mm	Fluted	Downward, Full	
3S45FUD	3mm	10mm	45°	12mm	Diamond	Downward, Partial	
3S45FI	3mm	10mm	45°	12mm	Fluted	Upward, Partial	
3L75FU	3mm	18mm	75°	19mm	Fluted	Downward, Partial	
2L75FU	2mm	18mm	75°	19mm	Fluted	Downward, Partial	
3L10DS	3mm	6mm	10°	12mm	Fluted	Distal, Full	



3S45FU



3S45FUS



3S45FUD



3S45FI



3L75FU



2L75FU



3L10DS

The Dreal™ is intended to cut bone in neurosurgical and spinal applications. Prior to use, the surgeon should refer to the product Instructions for Use (IFU) for complete indications, warnings, precautions and contraindications.

dreaTM Model selection guide

	Decompression (Above the root)	Decompression (Below the root)	TLIF/PLIF	ALIF	Cervical foraminotomy; uncusectomy	Osteophyte removal (ACDF)	Osteophyte removal (following corpectomy)	Osteophyte removal (lumbar)
Approach	Posterior / Lateral	Posterior / Lateral	Posterior / Lateral	Anterior	Anterior	Anterior	Anterior	Posterior
3S45FU	●			●	●			
3S45FUS	●							
3S45FUD	●			●	●			
3S45FI		●	●					●
3L75FU	●			●				
2L75FU							●	
3L10DS						●		

